


E_I-I: The Frequency of Vibration of a Stretched String: Variation with Length and Tension

Apparatus

Sonometer with two bridges; set of tuning forks (frequencies 256 to 512 Hz); 5m sample of sonometer wire; tuning fork hammer; sounding box; beam balance with masses; 2 sheets graph paper



Precaution: do not exceed the maximum tension (see label beside force meter). Place bridge A near the left hand end. Do not move it again.

Experiment 1: Procedure

To study how frequency varies with length and thus calculate the mass per unit length of the sonometer wire.

1. Set the tension to 20 N; throughout this experiment ensure that this remains constant. Record this value of tension.
2. Strike the 256 Hz tuning fork with the hammer, and press the tip hard against the sounding box (do not touch the arms of the fork). Listen carefully to the note produced.
3. Pluck the string between bridges A and B, then move bridge B until the string

makes exactly the same note as the fork. (Repeat step 2 to help the comparison).

4. Record the frequency and the length AB .
5. Repeat steps 2, 3 and 4, using all the other forks in turn, up to 512Hz.

Experiment 1: Theory

For a string under tension, secured at both ends:

$$f = \frac{1}{2l} \sqrt{\frac{T}{m}}$$

Where:

f = frequency (Hz)

m = mass per unit length (kg m^{-1})

l = length (m)

T = tension (N)

If T and m are constant, then $f = k \left(\frac{1}{l} \right)$, where the constant k is $\frac{1}{2} \sqrt{\frac{T}{m}}$.

Experiment 1: Analysis

1. Plot a graph of f against $\frac{1}{l}$ and find the gradient.
2. Use the gradient, the value of tension and the theory to find the value of m -- the mass per unit length of the wire.

Experiment 2: Procedure

To study how frequency varies with tension, and thus calculate the mass per unit length of the wire.

1. Place bridge B so that length AB is 20cm. Record this length.
2. Adjust the tension so that the string frequency is 256Hz (use the tuning fork as described in experiment 1 procedure step 2.).

3. Record the values of tension and frequency.
4. Adjust the tension so that the string produces the frequencies of each of the other tuning forks. For each fork record the values of tension and frequency. Ensure that length AB remains constant throughout.

NB: DO NOT EXCEED THE MAXIMUM TENSION ALLOWED

Experiment 2: Theory

In this experiment, l and m are constant, therefore:

$$f = k\sqrt{T} \quad \text{where the constant } k = \frac{1}{2l\sqrt{m}}$$

Experiment 2: Analysis

1. Plot a graph of f against \sqrt{T} , and find the gradient.
2. Use the gradient, the value of length, and the theory to find the value of m -- the mass per unit length of the wire.

Conclusion

Use the beam balance to measure the mass of the 5m length of sonometer wire supplied. From this value calculate the mass per unit length of the wire. Assuming that this value is very accurate, calculate the % error in the values obtained from experiments 1 and 2.

© 2015 [CC-BY](#) by Bob Drach and Norman Price
Based off of book published ????

[About](#)